


Anti-Gpkow Rabbit Polyclonal Antibody For Wb And Elisa - Q92917

Item Number: CM0015150



Introduction

Rabbit polyclonal antibody targeting human GPKOW (G-patch domain and KOW motifs-containing protein). Validated for WB and ELISA; cross-reacts with mouse, rat. Involved in pre-mRNA splicing. - Q92917

[Learn More](#)

Parameter	Value
Product Name	GPKOW Rabbit pAb
Remarks/Alias	T54; Mos2; Spp2; GPATC5; GPATCH5; GPKOW
Species	Human
Gene ID (Human)	27238
Gene ID	27238
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-476 of human GPKOW (NP_056513.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92917
Protein Weight	52kDa
Shipping	Ice bag