

Anti-Gpld1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P80108

Item Number: CM0025351

Introduction

GPLD1 Rabbit pAb validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat GPI-specific phospholipase D. Ideal for GPI anchor metabolism and signaling research.

[Learn More](#)

Parameter	Value
Product Name	GPLD1 Rabbit pAb
Remarks/Alias	PLD, GPIPLD, PIGPLD, GPIPLDM, PIGPLD1, GPLD1
Species	Human
GeneID (Human)	2822
GeneID	2822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-160 of human GPLD1 (NP_001494.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P80108
Protein Weight	92kDa
Shipping	Ice bag