

Anti-Gpm6A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P51674

Item Number: CM0027902

Introduction

Rabbit polyclonal antibody targeting human GPM6A (M6A), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. For neuronal differentiation and plasticity research.

[Learn More](#)

Parameter	Value
Product Name	GPM6A Rabbit pAb
Remarks/Alias	M6A; GPM6; GPM6A
Species	Human
GeneID (Human)	2823
GeneID	2823
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPM6A (NP_005268.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51674
Protein Weight	31kDa
Shipping	Ice bag