

Anti-Gpn1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcn4

Item Number: CM0014243

Introduction

High-quality rabbit polyclonal antibody targeting human GPN1 (XAB1/MBDIN). Validated for WB and ELISA, cross-reacts with mouse. Ideal for RNA polymerase II assembly and nuclear import studies.

[Learn More](#)

Parameter	Value
Product Name	GPN1 Rabbit pAb
Remarks/Alias	XAB1; MBDIN; NTPBP; RPAP4; ATPBD1A; GPN1
Species	Human
Gene ID	11321
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-374 of human GPN1 (NP_009197.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9HCN4
Protein Weight	42kDa
Shipping	Ice bag