

Anti-Gpr108 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npr9

Item Number: CM0020990

Introduction

CamelBio offers anti-GPR108 rabbit polyclonal antibody (LUSTR2) for WB and ELISA. Validated in human, mouse, and rat samples. Suitable for Golgi and AAV transduction research. SWISS: Q9NPR9.

[Learn More](#)

Parameter	Value
Product Name	GPR108 Rabbit pAb
Remarks/Alias	LUSTR2; GPR108
Species	Human
GeneID (Human)	56927
GeneID	56927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-260 of human GPR108 (NP_001073921.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPR9
Protein Weight	61kDa
Shipping	Ice bag