

Anti-Gpr109A/Hm74A/Hcar2 Monoclonal Antibody For Wb, Elisa - Q8Tds4

Item Number: CM0004041



Introduction

Rabbit anti-GPR109A/HCAR2 monoclonal antibody validated for WB and ELISA; cross-reacts with mouse and rat. Targets the high affinity niacin receptor that regulates adiponectin secretion and lipolysis, suitable for metabolic and inflammatory research.

[Learn More](#)

Parameter	Value
Product Name	GPR109A/HM74A/HCAR2 Rabbit mAb
Remarks/Alias	HCA2; HM74a; HM74b; PUMAG; NIACR1; Puma-g; GPR109A; GPR109A/HM74A/HCAR2
Species	Human
GeneID (Human)	338442
GeneID	338442
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human GPR109A/HM74A/HCAR2 (NP_808219.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8TDS4
Protein Weight	42kDa
Shipping	Ice bag