

Anti-Gpr141 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z602

Item Number: CM0016596

Introduction

Rabbit polyclonal antibody against human GPR141, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for orphan GPCR research and cell membrane receptor studies.

[Learn More](#)

Parameter	Value
Product Name	GPR141 Rabbit pAb
Remarks/Alias	PGR13; GPR141
Species	Human
GeneID	353345
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 50-150 of human GPR141 (NP_001316922.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7Z602
Protein Weight	35kDa
Shipping	Ice bag