

Anti-Gpr151 Polyclonal Antibody For Wb, Elisa - Q8Tdv0

Item Number: CM0020851

Introduction

Rabbit polyclonal antibody targeting human GPR151, a proton-sensing GPCR. Validated for Western blot and ELISA with cross-reactivity to mouse and rat. Ideal for neuroscience and GPCR studies.

[Learn More](#)

Parameter	Value
Product Name	GPR151 Rabbit pAb
Remarks/Alias	GPCR; PGR7; GALR4; GALRL; GPCR-2037; GPR151
Species	Human
GeneID (human)	134391
GeneID	134391
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-419 of human GPR151 (NP_919227.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDV0
Protein Weight	47kDa
Shipping	Ice bag