


Anti-Gpr155 Polyclonal Antibody For Wb, Elisa - Q7Z3F1

Item Number: CM0024507



Introduction

Rabbit polyclonal antibody targeting human GPR155, a lysosomal cholesterol sensor for mTORC1 signaling. Validated for Western blot and ELISA, cross-reactive with mouse and rat. Suitable for metabolic signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	GPR155 Rabbit pAb
Remarks/Alias	DEP.7; PGR22; DEPDC3; GPR155
Species	Human
GeneID (Human)	151556
GeneID	151556
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 721-870 of human GPR155 (NP_001028217.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7Z3F1
Protein Weight	97kDa
Shipping	Ice bag