


Anti-Gpr17 Rabbit Polyclonal Antibody For Wb, Elisa - Q13304

Item Number: CM0030332



Introduction

CamelBio offers GPR17 rabbit polyclonal antibody (Q13304) validated for WB, ELISA, detecting human, mouse, rat. Ideal for uracil nucleotide/cysteinyl leukotriene receptor research and brain ischemia studies.

[Learn More](#)

Parameter	Value
Product Name	GPR17 Rabbit pAb
Remarks/Alias	GPR17
Species	Human
Gene ID	2840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human GPR17 (NP_005282.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13304
Protein Weight	41kDa
Shipping	Ice bag