

Anti-Gpr173 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ns66

Item Number: CM0016155

Introduction

Rabbit polyclonal antibody against GPR173 (SREB3), validated for WB and ELISA in human and mouse samples. Targets the phoenixin receptor involved in GNRH signaling and granulosa cell proliferation; ideal for neuroendocrine and reproductive research.

[Learn More](#)

Parameter	Value
Product Name	GPR173 Rabbit pAb
Remarks/Alias	SREB3; GPR173
Species	Human
GeneID (Human)	54328
GeneID	54328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-230 of human GPR173 (NP_061842.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NS66
Protein Weight	41kDa
Shipping	Ice bag