

Anti-Gpr22 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99680

Item Number: CM0022147

Introduction

GPR22 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat GPR22. Orphan GPCR involved in G(i)/G(o) signaling and ciliogenesis. Ideal for cardiovascular and neuroscience research applications.

[Learn More](#)

Parameter	Value
Product Name	GPR22 Rabbit pAb
Remarks/Alias	GPR22
Species	Human
GeneID (Human)	2845
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPR22 (NP_005286.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99680
Protein Weight	49kDa
Shipping	Ice bag