

Anti-Gpr32 Polyclonal Antibody For Wb, Elisa - 075388

Item Number: CM0023288

Introduction

Rabbit polyclonal antibody against human GPR32 (DRV1/RVDR1) for Western blot and ELISA. Validated cross-reactivity with human, mouse, and rat. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	GPR32 Rabbit pAb
Remarks/Alias	DRV1; RVDR1; GPR32
Species	Human
GeneID (Human)	2854
GeneID	2854
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPR32 (NP_001497.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75388
Protein Weight	40kDa
Shipping	Ice bag