

Anti-Gpr35 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hc97

Item Number: CM0024117



Introduction

Rabbit polyclonal antibody targeting human GPR35, a GPCR involved in immune response, cardioprotection, and inflammation. Validated for Western blot and ELISA, with cross-reactivity to mouse. Ideal for GPCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	GPR35 Rabbit pAb
Remarks/Alias	GPR35
Species	Human
Gene ID (Human)	2859
Gene ID	2859
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human GPR35 (NP_005292.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9HC97
Protein Weight	34kDa
Shipping	Ice bag