

Anti-Gpr37 Rabbit Polyclonal Antibody For Wb And Elisa - O15354

Item Number: CM0022038

Introduction

Rabbit polyclonal antibody targeting the prosaposin receptor GPR37 (PAELR/EDNRBL), suitable for Western blot and ELISA. Detects human and rat GPR37, supporting studies of myelination, inflammation, and neurodegeneration.

[Learn More](#)

Parameter	Description
Product Name	GPR37 Rabbit pAb
Remarks/Alias	PAELR; EDNRBL; hET(B)R-LP; GPR37
Species	Human
GeneID (human)	2861
GeneID	2861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-180 of human GPR37 (NP_005293.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O15354
Protein Weight	67kDa
Shipping	Ice bag