

Anti-Gpr4 Polyclonal Antibody For Wb, Elisa - P46093

Item Number: CM0021041



Introduction

High-quality rabbit polyclonal antibody against human GPR4, validated for WB and ELISA. Suitable for studying proton-sensing GPCR signaling in pH homeostasis and respiratory control.

[Learn More](#)

Parameter	Value
Product Name	GPR4 Rabbit pAb
Remarks/Alias	GPR6C.I; GPR4
Species	Human
GeneID (Human)	2828
GeneID	2828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 283-362 of human GPR4 (NP_005273.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P46093
Protein Weight	41kDa
Shipping	Ice bag