

Anti-Gpr55 Polyclonal Antibody For Wb, If-P, Elisa - Q9Y2T6

Item Number: CM0016482



Introduction

High-quality rabbit polyclonal antibody against human GPR55 (LPIR1), validated in WB, IF-P, ELISA, with cross-reactivity to mouse and rat. For research use in GPCR signaling and lipid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	GPR55 Rabbit pAb
Remarks/Alias	LPIR1; GPR55
Species	Human
GeneID (Human)	9290
GeneID	9290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-319 of human GPR55 (NP_005674.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y2T6
Protein Weight	37kDa
Shipping	Ice bag