


Anti-Gpr6 Polyclonal Antibody For Wb, Elisa - P46095

Item Number: CM0017378



Introduction

High-quality rabbit polyclonal antibody targeting human GPR6 for Western blot and ELISA. Detects endogenous GPR6 in human and mouse samples. Suitable for GPCR signaling and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	GPR6 Rabbit pAb
Remarks/Alias	GPR6
Species	Human
Gene ID (Human)	2830
Gene ID	2830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 253-362 of human GPR6 (NP_005275.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46095
Protein Weight	38kDa
Shipping	Ice bag