

Anti-Gpr62 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bzj7

Item Number: CM0024513

Introduction

Rabbit polyclonal antibody targeting human GPR62 (GPCR8, KPG_005), raised against a synthetic peptide. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for studying orphan GPCR signaling in brain tissues.

[Learn More](#)

Parameter	Specification
Product Name	GPR62 Rabbit pAb
Remarks/Alias	GPCR8; KPG_005; GPR62
Species	Human
GeneID	118442
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human GPR62 (NP_543141.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZJ7
Protein Weight	38kDa
Shipping	Ice bag