


Anti-Gpr68 Polyclonal Antibody For Wb, Elisa - Q15743

Item Number: CM0019388



Introduction

Rabbit polyclonal antibody against human GPR68, validated for WB and ELISA. Detects endogenous GPR68 in human, mouse samples. Ideal for pH sensing and mechanotransduction research.

[Learn More](#)

Parameter	Value
Product Name	GPR68 Rabbit pAb
Remarks/Alias	OGR1; AI2A6; GPR12A; GPR68
Species	Human
GeneID (Human)	8111
GeneID	8111
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 276-365 of human GPR68 (NP_003476.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15743
Protein Weight	41kDa
Shipping	Ice bag