

Anti-Gpr82 Rabbit Polyclonal Antibody For Wb, Elisa - Q96P67

Item Number: CM0019111

Introduction

Investigate the orphan G protein-coupled receptor 82 with this rabbit polyclonal antibody. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cell membrane receptor studies.

[Learn More](#)

Parameter	Value
Product Name	GPR82 Rabbit pAb
Remarks/Alias	GPR82
Species	Human
GeneID	27197 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPR82 (NP_543007.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96P67
Protein Weight	38kDa
Shipping	Ice bag