

Anti-Gprc5A Polyclonal Antibody For Wb, Elisa - Q8Nfj5

Item Number: CM0027247

Introduction

Rabbit polyclonal antibody targeting GPRC5A (RAIG-1), validated for Western blot and ELISA. Recognizes human and mouse retinoic acid-induced protein 3, a lung tumor suppressor. Ideal for cancer research and GPCR studies.

[Learn More](#)

Parameter	Value
Product Name	GPRC5A Rabbit pAb
Remarks/Alias	RAI3; TIG1; RAIG1; GPCR5A; PEIG-1; GPRC5A
Species	Human
Gene ID (Human)	9052
Gene ID	9052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 268-357 of human GPRC5A (NP_003970.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NFJ5
Protein Weight	40kDa
Shipping	Ice bag