


Anti-Gprc5D Monoclonal Antibody For Fc - Q9Nzd1

Item Number: CM0012789



Introduction

Research-grade rabbit anti-human GPRC5D monoclonal antibody, -conjugated, for flow cytometry (FC). Validated on human samples; shipped on ice bags. UniProt: Q9NZD1. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GPRC5D mAb
Remarks/Alias	GPRC5D
Species	Human
Gene ID (Human)	55507
Gene ID	55507
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPRC5D (NP_061124.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NZD1
Protein Weight	34kDa/36kDa/38kDa
Shipping	Ice bag