

Anti-Gprc6A Rabbit Polyclonal Antibody For Wb, Elisa - Q5T6X5

Item Number: CM0026463

Introduction

Rabbit polyclonal antibody targeting human GPRC6A (GPCR33) for WB and ELISA. Cross-reacts with mouse and rat. Suitable for investigating GPRC6A-mediated signaling, hormone regulation, and tissue expression.

[Learn More](#)

Parameter	Value
Product Name	GPRC6A Rabbit pAb
Remarks/Alias	GPCR; bA86F4.3; GPRC6A
Species	Human
GeneID (Human)	222545
GeneID	222545
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-590 of human GPRC6A (NP_683766.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T6X5
Protein Weight	105kDa
Shipping	Ice bag