


Anti-Gpsm2 Rabbit Polyclonal Antibody For Wb, Elisa - P81274

Item Number: CM0013761



Introduction

High-quality GPSM2 rabbit polyclonal antibody validated for WB and ELISA applications. Reacts with human and mouse samples; targets G protein-signaling modulator 2 (LGN) involved in mitotic spindle organization and asymmetric cell division. Research use only.

[Learn More](#)

Parameter	Details
Product Name	GPSM2 Rabbit pAb
Remarks/Alias	LGN; CMCS; PINS; DFN82; GPSM2
Species	Human
GeneID (Human)	29899
GeneID	29899
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-490 of human GPSM2 (NP_037428.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P81274
Protein Weight	77kDa
Shipping	Ice bag