

Anti-Gpt Rabbit Polyclonal Antibody For Wb, Elisa - P24298

Item Number: CM0024822

Introduction

Rabbit polyclonal antibody targeting human alanine aminotransferase 1 (GPT/ALT1). Suitable for Western blot and ELISA applications, with cross-reactivity to mouse and rat GPT. Useful for hepatic metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GPT Rabbit pAb
Remarks/Alias	ALT; AAT1; ALT1; GPT1; SGPT; GPT
Species	Human
GeneID (Human)	2875
GeneID	2875
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human GPT (NP_005300.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24298
Protein Weight	55kDa
Shipping	Ice bag