

Anti-Gpt2 Mouse Monoclonal Antibody For Wb - Q8Td30

Item Number: CM0007258



Introduction

Mouse monoclonal antibody against human Alanine aminotransferase 2 (GPT2/ALT2), raised against recombinant protein (aa 358-456). Validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Suitable for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GPT2 Mouse mAb
Remarks/Alias	ALT2; GPT 2; MRT49; NEDSPM; GPT2
Species	Human
GeneID (Human)	84706
GeneID	84706
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence within amino acids 358 - 456 of human GPT2 (Q8TD30).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TD30
Protein Weight	58kDa
Shipping	Ice bag