

Anti-Gpt2 Polyclonal Antibody For Wb, Elisa - Q8Td30

Item Number: CM0021895



Introduction

Rabbit anti-GPT2 polyclonal antibody validated for WB and ELISA. Reacts specifically with human GPT2 (alanine aminotransferase 2), a key metabolic enzyme with implications in neurological disorders. Ideal for amino acid metabolism research and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	GPT2 Rabbit pAb
Remarks/Alias	ALT2; GPT 2; MRT49; NEDSPM; GPT2
Species	Human
GeneID (Human)	84706
GeneID	84706
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 310-410 of human GPT2. (NP_597700.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TD30
Protein Weight	58kDa
Shipping	Ice bag