

Anti-Gpvi/Gp6 Rabbit Polyclonal Antibody For Wb, Elisa - P0C191

Item Number: CM0023914

Introduction

Rabbit polyclonal antibody targeting mouse GPVI/GP6, a platelet glycoprotein receptor for collagen. Validated for WB and ELISA with human cross-reactivity. Useful for studying platelet activation and thrombosis.

[Learn More](#)

Parameter	Value
Product Name	GPVI/GP6 Rabbit pAb
Remarks/Alias	Gpvi; Gm469; 9830166G18Rik; GPVI/GP6
Species	Mouse
GeneID	243816
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-265 of mouse GPVI/GP6 (NP_001156486.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P0C191
Protein Weight	34kDa
Shipping	Ice bag