


Anti-Gpx1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P07203

Item Number: CM0006372



Introduction

Rabbit monoclonal antibody against human GPX1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for studying redox homeostasis and hydroperoxide reduction.

[Learn More](#)

Parameter	Value
Product Name	GPX1 Rabbit mAb
Remarks/Alias	GPXD; GSHPX1; GPX1
Species	Human
GeneID (Human)	2876
GeneID	2876
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human GPX1 (NP_000572.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07203
Protein Weight	22kDa
Shipping	Ice bag