


Anti-Gpx1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P07203

Item Number: CM0014697



Introduction

Rabbit polyclonal antibody against human GPX1 for WB, IF/ICC, IHC-P, ELISA. Recognizes human, mouse, rat GPX1. Ideal for redox homeostasis and lipid peroxidation studies.

[Learn More](#)

Parameter	Value
Product Name	GPX1 Rabbit pAb
Remarks/Alias	GPXD, GSHPX1, GPX1
Species	Human
GeneID	2876
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human GPX1 (NP_000572.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07203
Protein Weight	22 kDa
Shipping	Ice bag