


Anti-Gpx2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P18283

Item Number: CM0014193



Introduction

GPX2 Rabbit pAb detects human, mouse, and rat GPX2 (P18283) via WB, IHC-P, ELISA. Recombinant immunogen, validated for redox and GI research applications. Optimized for studies of glutathione peroxidase 2 function in liver and gastrointestinal tissues, this polyclonal antibody ensures reliable detection in multiple assay formats.

[Learn More](#)

Parameter	Value
Product Name	GPX2 Rabbit pAb
Remarks/Alias	GPRP; GPx-2; GI-GPx; GPRP-2; GPx-GI; GSHPx-2; GSHPX-GI; GPX2
Species	Human
GeneID (Human)	2877
GeneID	2877
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human GPX2 (NP_002074.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18283
Protein Weight	22kDa
Shipping	Ice bag