


Anti-Gpx3 Rabbit Polyclonal Antibody For Wb And Elisa - P22352

Item Number: CM0018475



Introduction

Human GPX3 rabbit polyclonal antibody (GPx-P) for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Targets extracellular glutathione peroxidase, an antioxidant enzyme protecting cells from oxidative damage.

[Learn More](#)

Parameter	Value
Product Name	GPX3 Rabbit pAb
Remarks/Alias	GPx-P; GSHPx-3; GSHPx-P; GPX3
Species	Human
GeneID (Human)	2878
GeneID	2878
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPX3 (NP_002075.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P22352
Protein Weight	26kDa
Shipping	Ice bag