


Anti-Gpx4 Monoclonal Antibody For Wb - P36969

Item Number: CM0002193



Introduction

Validate GPX4 protein expression with this knockout-validated rabbit monoclonal antibody. Suitable for WB, IF-P, ELISA applications across human, mouse, and rat samples. Ideal for ferroptosis and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GPX4 Rabbit mAb
Remarks/Alias	MCSP; SMDS; GPx-4; PHGPx; snGPx; GSHPx-4; snPHGPx; [KD Validated] GPX4
Species	Human
Gene ID (Human)	2879
Gene ID	2879
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPX4 (P36969).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36969
Protein Weight	22kDa
Shipping	Ice bag