


Anti-Gpx4 Rabbit Monoclonal Antibody (Kd Validated) For Wb, If-P, Ihc-P, Elisa - P36969

Item Number: CM0009232



Introduction

High-affinity anti-GPX4 rabbit monoclonal antibody (KD validated) for Western blot, IF, IHC, and ELISA. Cross-reacts with human, mouse, and rat GPX4. Ideal for ferroptosis and oxidative stress research.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] GPX4 Rabbit mAb
Remarks/Alias	MCSP; SMDS; GPx-4; PHGPx; snGPx; GSHPx-4; snPHGPx; [KD Validated] GPX4
Species	Human
Gene ID	2879
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 74-197 of human GPX4 (NP_002076.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36969
Protein Weight	22kDa
Shipping	Ice bag