


Anti-Gpx4 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ihc, If, Elisa - P36969

Item Number: CM0010583



Introduction

High-quality rabbit monoclonal antibody against human GPX4, validated for Western blot, IHC, IF, and ELISA. Detects a key antioxidant and ferroptosis regulator, with cross-reactivity across human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] GPX4 Rabbit mAb
Remarks/Alias	MCSP; SMDS; GPx-4; PHGPx; snGPx; GSHPx-4; snPHGPx
Species	Human
GeneID (Human)	2879
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 74-197 of human GPX4 (NP_002076.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36969
Protein Weight	20 kDa/22 kDa
Shipping	Ice bag