


Anti-Gpx4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36969

Item Number: CM0013042



Introduction

GPX4 Rabbit Polyclonal Antibody validates in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat GPX4, an essential phospholipid hydroperoxidase protecting cells from ferroptosis and oxidative damage. Ideal for redox biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	GPX4 Rabbit pAb
Remarks/Alias	MCSP; SMDS; GPx-4; PHGPx; snGPx; GSHPx-4; snPHGPx
Species	Human
GeneID (Human)	2879
GeneID	2879
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-197 of GPX4 (NP_002076.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36969
Protein Weight	22 kDa
Shipping	Ice bag