


Anti-Gpx7 Rabbit Polyclonal Antibody For Wb, Elisa - Q96SI4

Item Number: CM0025979



Introduction

Anti-GPX7 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat GPX7, an ER thiol peroxidase involved in redox homeostasis and oxidative protein folding. Ideal for oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	GPX7 Rabbit pAb
Remarks/Alias	GPX6; CL683; GPx-7; NPGPx; GSHPx-7; GPX7
Species	Human
GeneID (Human)	2882
GeneID	2882
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 20-187 of human GPX7 (NP_056511.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SL4
Protein Weight	21kDa
Shipping	Ice bag