

Anti-Gpx8 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ted1

Item Number: CM0022315



Introduction

Anti-GPX8 rabbit polyclonal antibody for WB and ELISA applications. Recognizes human, mouse, and rat GPX8, an ER thiol peroxidase crucial for redox homeostasis and disulfide formation. Supports oxidative stress studies.

[Learn More](#)

Parameter	Value
Product Name	GPX8 Rabbit pAb
Remarks/Alias	GPx-8; UNQ847; EPLA847; GSHPx-8; GPX8
Species	Human
Gene ID (Human)	493869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-209 of human GPX8 (NP_001008398.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TED1
Protein Weight	24kDa
Shipping	Ice bag