

Anti-Granzyme B Polyclonal Antibody For Wb, Elisa - P04187

Item Number: CM0028626

Introduction

Rabbit polyclonal antibody targeting mouse Granzyme B (P04187), suitable for WB and ELISA. Granzyme B is a key serine protease in cytotoxic T-cell and NK-cell-mediated apoptosis and pyroptosis.

[Learn More](#)

Parameter	Value
Product Name	Granzyme B Rabbit pAb
Remarks/Alias	GZB; CCP1; Ctla1; Ctla-1; CCP-1/C11
Species	Mouse
GeneID (human)	3002
GeneID	14939
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 141-240 of mouse Granzyme B (NP_038570.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04187
Protein Weight	27 kDa
Shipping	Ice bag