

Anti-Granzyme B Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P10144

Item Number: CM0023771



Introduction

Granzyme B rabbit polyclonal antibody for WB, IF-P, IHC-P, ELISA. Detects human and mouse Granzyme B (~28kDa). Ideal for apoptosis, pyroptosis, and immune cytotoxicity research. Synthetic peptide immunogen.

[Learn More](#)

Parameter	Value
Product Name	Granzyme B Rabbit pAb
Remarks/Alias	C11; HLP; CCPI; CGL1; CSPB; SECT; CGL-1; CSP-B; CTLA1; CTSG1; Granzyme B
Species	Human
Gene ID (Human)	3002
Gene ID	3002
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10144
Protein Weight	28kDa
Shipping	Ice bag