

Anti-Grap Rabbit Polyclonal Antibody For Wb And Elisa - Q13588

Item Number: CM0014804

Introduction

High-quality rabbit polyclonal antibody targeting human GRAP (GRB2-related adapter protein), validated for Western blot and ELISA. Reliable tool for Ras signaling, inner ear biology, and deafness research.

[Learn More](#)

Parameter	Specification
Product Name	GRAP Rabbit pAb
Remarks/Alias	DFNB114; GRAP
Species	Human
GeneID (Human)	10750
GeneID	10750
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 128-217 of human GRAP (NP_006604.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13588
Protein Weight	25kDa
Shipping	Ice bag