

Anti-Grasp65 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Bqq3

Item Number: CM0006630



Introduction

High-quality GRASP65 rabbit monoclonal antibody for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for Golgi biology and cell structure research.

[Learn More](#)

Parameter	Value
Product Name	GRASP65 Rabbit mAb
Remarks/Alias	P65; GOLPH5; GRASP65
Species	Human
GeneID	64689
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQQ3
Protein Weight	46kDa
Shipping	Ice bag