

Anti-Grasp65 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bqq3

Item Number: CM0029243



Introduction

Rabbit polyclonal antibody against human GRASP65 (GORASP1, GOLPH5), validated for WB, IF/ICC, and ELISA. Detects a 46 kDa Golgi stacking protein involved in Golgi ribbon formation, mitosis, and vesicle docking. Ideal for apoptosis and trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	GRASP65 Rabbit pAb
Remarks/Alias	P65; GOLPH5; GRASP65
Species	Human
Gene ID (Human)	64689
Gene ID	64689
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 221-440 of human GRASP65 (NP_114105.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9BQQ3
Protein Weight	46kDa
Shipping	Ice bag