


Anti-Grb2 Polyclonal Antibody For Wb, If/Icc, Elisa - P62993

Item Number: CM0025884



Introduction

High-quality anti-GRB2 polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human and mouse GRB2, a key adapter in signaling pathways. Ideal for studying cell signaling, immunity, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	GRB2 Rabbit pAb
Remarks/Alias	ASH; Grb3-3; MST084; NCKAP2; MSTP084; EGFRBP-GRB2; GRB2
Species	Human
GeneID (Human)	2885
GeneID	2885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-217 of human GRB2 (NP_002077.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P62993
Protein Weight	25kDa
Shipping	Ice bag