

Anti-Grb7 Rabbit Polyclonal Antibody For Wb And Elisa - Q14451

Item Number: CM0029426



Introduction

Rabbit polyclonal antibody targeting human GRB7, validated for WB and ELISA. Suitable for investigating GRB7-mediated signaling, cell proliferation, migration, and stress granule dynamics in research and IVD applications.

[Learn More](#)

Parameter	Value
Product Name	GRB7 Rabbit pAb
Remarks/Alias	GRB7
Species	Human
Gene ID (Human)	2886
Gene ID	2886
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human GRB7 (NP_001229372.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14451
Protein Weight	60kDa
Shipping	Ice bag