


Anti-Grem1 Polyclonal Antibody For Wb, Elisa - O60565

Item Number: CM0015247



Introduction

High-quality rabbit polyclonal antibody against human GREM1 (Gremlin-1), validated for Western blot and ELISA. Ideal for research on BMP antagonism, carcinogenesis, and kidney organogenesis.

[Learn More](#)

Parameter	Value
Product Name	GREM1 Rabbit pAb
Remarks/Alias	DRM; HMPS; MPSH; PIG2; CRAC1; CRCS4; DAND2; HMPS1; IHG-2; DUP15q; C15DUPq; GREMLIN; CKTSF1B1; GREM1
Species	Human
GeneID (Human)	26585
GeneID	26585
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-143 of human GREM1 (NP_001178252.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60565
Protein Weight	21kDa
Shipping	Ice bag