

Anti-Grhl2 Polyclonal Antibody For Wb And Elisa - Q6Isb3

Item Number: CM0018029



Introduction

Rabbit anti-GRHL2 polyclonal antibody validated for WB and ELISA; reacts with human and rat. Targets transcription factor GRHL2 (71 kDa), critical for epithelial development, neurulation, and tight junction regulation. Ideal for developmental and epithelial research.

[Learn More](#)

Parameter	Value
Product Name	GRHL2 Rabbit pAb
Remarks/Alias	BOM, ECTDS, PPCD4, DFNA28, TFCP2L3, GRHL2
Species	Human
GeneID (Human)	79977
GeneID	79977
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-250 of human GRHL2 (NP_079191.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6ISB3
Protein Weight	71kDa
Shipping	Ice bag