

Anti-Grhpr Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubq7

Item Number: CM0024000



Introduction

Rabbit polyclonal antibody targeting GRHPR, validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying glyoxylate metabolism and hydroxypyruvate reductase enzyme activity.

[Learn More](#)

Parameter	Value
Product Name	GRHPR Rabbit pAb
Remarks/Alias	PH2; GLXR; GLYD; GRHPR
Species	Human
GeneID (Human)	9380
GeneID	9380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-328 of human GRHPR (NP_036335.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBQ7
Protein Weight	36kDa
Shipping	Ice bag