

Anti-Gria3 Polyclonal Antibody For Wb, If/Icc, Elisa - P42263

Item Number: CM0016209

Introduction

GRIA3 Rabbit pAb is a polyclonal antibody raised in rabbit against human GRIA3, validated for WB, IF/ICC, ELISA. It detects human, mouse, rat GRIA3, a 101 kDa AMPA-selective glutamate receptor involved in excitatory synaptic transmission. Ideal for neuroscience research.

[Learn More](#)

Parameter	Description
Product Name	GRIA3 Rabbit pAb
Remarks/Alias	GLUR3; GLURC; GluA3; MRX94; MRXSW; GLUR-C; iGluR3; GLUR-K3; GRIA3
Species	Human
GeneID	2892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-490 of human GRIA3 (NP_000819.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42263
Protein Weight	101kDa
Shipping	Ice bag